Weropress®

Kaltpolymerisat

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Weropress®

Cold-curing polymer

(en) Instructions for use, please read carefully!

(USA) Caution: Federal law restricts this device to sale by or on the order of a licensed dental practitioner.

Cold-curing, colour-stable, methylmethacrylate-based resin for fabricating partial and total prostheses using injection, tamping-pressing and casting processes.

Composition

Weropress® powder contains

- PMMA (polymer and copolymer)
- barbituric acid catalyst system
- organic colorants
- inorganic pigments

Weropress® liquid contains

- MMA (methylmethacrylate)
- dimethacrylate
- barbituric acid catalyst system

Dosage and application

- See table for mixing ratios and processing times depending on the processing method of the acrylic.
- Always wet casts sufficiently for 5 to 15 minutes depending on the degree of dryness.
- Use standard alginate isolation to separate plaster from acrylic.
- Do not compress the powder when using dosage aids.

Injection process

- Proceed as specified by the device manufacturer
- Use class IV plaster for making casts and investment where possible.
- The temperature of the plaster in the flask should be 35 - 40 °C / 95 - 104 °F to promote polymerisation.

Press process


- Use class IV plaster for making casts and investment where possible.
- The temperature of the plaster in the flask should be 35 - 40 °C / 95 - 104 °F to promote polymerisation.
- Complete the press process promptly under the hydraulic press within 1 minute.

Casting process

- Use class IV plaster for making casts where possible.
- Investment with plaster, silicone or hydrocolloid (gel) depending on the process.
- Polymerisation is conducted at 2 to 2.5 bar pressure in water at 45 °C / 113 °F after appropriate dwell time.

Times and quantities

	Injection process	Tamping-pressing process	Casting process
Mixing ratio*	30 g : 15 ml	30 g : 15 ml	30 g : 21 ml
Polymerisation temperature	Residual heat in the flask	Residual heat in the flask	45 °C / 113 °F
Polymerisation time	15 min	15 min	15 min ¹ / 30 min ²
Pressure required during polymerisation	By device type	Pressure in the flask must be maintained during polymerisation	2 - 2,5 bar

Processing the acrylic

To reduce the residual monomer content, the finished (sanded) but not polished polymerised workpiece should be stored in water at room temperature for about 48 hours before delivery.

Trimming

To prevent inaccurate fitting after polymerisation, avoid excessive heat development during trimming and polishing.

Bond

- Roughen the basal surfaces of the acrylic teeth before bonding.
- Ensure that no impurities such as wax or isolating agent are present.
- Mechanical retention is recommended.
- Adequate mechanical retention is essential for ceramic teeth.

Restoration

Products of Weropress® can be restored and supplemented at any time with all standard MMA-based cold-curing polymers, such as Weropress®.

Notes

- Close all powder and liquid containers tightly after use.
- Store in a dry and cool area away from direct sunlight.
- Take precautionary measures against static discharges.
- Do not use the material after the expiration date.
- Record the LOT number with every process that requires identification of the material.
- The residual monomer content after polymerisation of Weropress® is < 2%.

Disposal information

Disposal of contents/container in accordance with local/regional/national/international regulations.

Danger and safety instruction Weropress® Monomer

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P272 Contaminated work clothing should not be allowed out of the workplace.

P262 Do not get in eyes, on skin, or on clothing.

P210 Keep away from heat. No Smoking.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

Danger and safety instructions Weropress® Polymer

H317 May cause an allergic skin reaction.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

Classification

Weropress® conforms to EN ISO 20795-1 Type 2 Class 1.

The product properties are based on compliance with the instructions for use.

For detailed information see the applicable Material Safety Data Sheet.

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Merz Dental GmbH