EC [REP]	SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, California 92078-4059, USA
SHOFU DENTAL GmbH Am Brüll 17, 40878 Ratingen, Germany	

SHOFU DENTAL ASIA-PACIFIC PTE LTD. 10 Science Park Road, #03-12, The Alpha, Science Park II, Singapore 117684	EN Instructions for use NL Gebruiksaanwijzing
--	--

DE Gebrauchsanleitung ES Instrucciones de uso	FR Mode d'emploi IT Istruzioni per l'uso	03:2019-03 71600
--	---	---------------------

EN	READ CAREFULLY BEFORE USE
-----------	----------------------------------

Light-curing C+B Composite

PRODUCT DESCRIPTION

SOLIDEX is an indirect veneering composite for metal based crowns and bridges, inlays/onlays, veneers, anterior jacket crowns as well as for temporary crowns and repair of indirect restorations intraorally and extraorally. Application of METAL PHOTO PRIMER or PRIMER PASTE to the metal surface enhances the bond strength between metal and SOLIDEX materials. For characterization or modification of the restorations, light-curing coloring paste, LITE ART, is recommended.

INDICATIONS FOR USE

➤ Metal based crowns and bridges
➤ Inlays/Onlays, Anterior jacket crowns, Laminate veneers
➤ Intraoral and extraoral repairs of metal based, composite or porcelain crowns

PRECAUTIONS AND WARNINGS

- Do not use this product on patients with known allergies to this material and/or methacrylate monomer.
- Operators with known allergies to this material and/or methacrylate monomer must not use this product.
- If any inflammation or other allergic reactions occur on either the patient operator, immediately discontinue use and seek medical advice.
- Do not this product on patients with malocclusion and bruxism (clenching, grinding and tapping).
- Do not handle uncured material of this product with bare hands. Use plastic gloves and protective glasses to avoid contact with skin and/or eyes. In case of accidental contact with skin, immediately blot with alcohol moistened cotton ball. In case of contact with eyes, immediately flush the eyes with plenty of water and seek medical advice.
- Uncured material of this product has caused inflammation or other allergic reactions when handled with bare hands.
- Restorations placed in patients’ mouths might stain or adhere plaque due to their dietary habits. Give the patient proper instructions for daily cleaning to avoid getting stained or plaque adhesion.
- Since SOLI BOND and METAL PHOTO PRIMER are inflammable, avoid any source of fire.
- Use this material in well ventilated areas (several appropriate ventilations per hour are recommended).
- Use local dust extractor, dust protective mask, etc while grinding this product to avoid harmful influence of the dust on the human body.
- Follow the instructions for use of each dental material, instrument or equipment used in conjunction with this product.
- Do not use this product for any purposes other than specifically outlined in these instructions for use.
- Use this product within the expiration date indicated on the package. (Example ☞ YYYY-MM-DD → Year-Month-Date of the expiration date)
- This product is intended for use by dental professionals only.

DIRECTIONS FOR USE

1. Metal framework

In order to obtain an optimal shade, it is important to design the metal framework according to the generally accepted guidelines. Apply retention beads (100-200 µm) on the wax pattern before casting. When using a silane bonding system, the application of retention beads can be dispensed in areas with limited space.

2. Metal treatment

Sandblast with 50-100 µm aluminum oxide and then clean with either steam or ultrasonic. Dry the surface with oil free air.

3. Bonding systems

Select (A) or (B) prior to the following opaque application.

(A) PRIMER PASTE - Primer paste is a metal bonding system and first layer of opaque in one product. The milky dull paste has to be applied evenly, that the retention beads are covered. Light cure for 1-3 minutes in a light curing unit (420-500 nm) tested and recommended by SHOFU INC.

(B) METAL PHOTO PRIMER - Apply a thin layer of metal primer over the surface to be veneered with composite. Then dry for 30 - 60 seconds.

(C) SOLI BOND- This is a silane coupling agent which is used in conjunction with silane bonding systems. e.g. Silicoater, Silicoater MD (Heraeus Kulzer), or Rocatec (3M ESPE).

4. Opaque

When using METAL PHOTO PRIMER, it is necessary to apply the first opaque (PRE-OPAQUE of SHOFU Universal Opaque) thinly. Light-cure for 3-5 minutes. Mask the framework completely with subsequent layers of shade opaque.

➤ When using SOLIDEX OPAQUE, apply a thin first layer, taking care that the retention beads are fully covered.

➤ Color modification of the opaque: If more translucency is required in the incisal or interproximal areas, InO or VrO can be used. Light-cure for 3-5 minutes.

5. Base Paste

For bridge units, it is recommended that BASE paste is applied after the first opaque layer to the base of the pontic area to the level of the adjacent crown. Light-cure for 3-5 minutes. Then, cover with the second opaque layer.

6. Cervical and Body

Apply Cervical and pre-cure for 1-3 minutes. Then apply Body into the correct anatomical form. Light-cure for 1-3 minutes.

7. Incisal and Translucent

Apply Incisal to complete the anatomical form, and light-cure for 3 minutes.

If Translucent is used, apply Incisal first and light-cure for 1-3 minutes. Then apply Translucent and finally light-cure for 3-5 minutes.

8. Adjusting

Polymerized surfaces must be roughened (e.g. sandblast with aluminum oxide before further pastes are applied). Scrub off the dust with a brush or compressed air. Do not clean with water or steam cleaner. Apply a thin coat of CERAMAGE MODELLING LIQUID.

9. Finishing and Polishing

Use silicon carbide burs (Dura-Green Stones) or diamond burs for contouring.

Then, polish with SoftCut PA/PB or AcryPoints PC2. Use DURA-POLISH with a brush for pre polishing and DURA-POLISH DIA with a leather/wool mops for high luster polishing.

NOTE ON USE

- When storing these materials in a refrigerator, allow them to reach room temperature before use.
- Tightly close the cap after each use.
- The paste dispensed onto the dish or paper pad should be shielded from ambient light by a light-shielding cover to avoid unintended polymerization of the material.
- Do not use this product in places with intense light such as close to a window or under the work light to avoid unintended polymerization of the paste.
- Base Paste is designed exclusively for pontics. Do not use for facing crowns.

COMPOSITION

PRIMER PASTE:	UDMA, 4-AET, Aluminum-Silicate, Glass powder and others
METAL PHOTO PRIMER:	Acetone, 4-AET and others
SOLI BOND:	Dehydrated ethanol, Silane coupling agent
OPAQUE:	UDMA, 4-AET, Aluminum Silicate, Organic filler Glass powder, Pigments and others
C&B COMPOSITE RESIN:	UDMA, Organic filler, Silica powder, Pigments and others

STORAGE

Store in a dark place at room temperature (1-30 °C/34-86 °F). Avoid high temperatures and high humidity. Keep away from direct sunlight and any source of ignition.

CAUTION: US Federal law restricts this device to sale by or on the order of a dental professional.

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