

Surface



Technical Data Sheet

Intended use: Dental sealing lacquer for coating resin temporaries.

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical term
1	Methyl methacrylate
2	Multifunctional urethane acrylate
3	Initiators
4	Additive

Properties:

Colour: clear transparent
Form: liquid
Density: 1.07 g cm³
Wave length polymerisation: 300 – 400 [nm]
Polymerisation time: 2-5 [min], depends on device
Storage: at room temperature, 18° - 28° [C],
keep the bottle closed

Clinical tests: Zahn-,Mund- und Kieferklinik WWU
(Westfälische Wilhelms-Universität Münster 1991)

Processing:

Apply a thin layer, provides a smooth and homogenous surface and increases the durability of the colouring inside the mouth.

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
dentamid@dreve.de
www.dreve.com